



100C • 100 • 150 • 250C • 250 • 250H • 350




USER MANUAL

CONTENTS

I. INTRODUCTION	4
II. DIRECTION FOR USE	6
1. Intended use	7
a. Intended purpose	7
2. Expected clinical benefit	7
3. Contraindications	7
4. Side effects	7
5. Intended population	7
6. Intended users	7
III. CAUTIONS AND WARNINGS	8
1. Definitions	9
2. Product safety	9
a. Precautions for use	9
b. Power source	10
c. Electromagnetic compatibility	10
d. Others	10
IV. PRODUCT DESCRIPTION	11
1. Product plan with description	12
a. OST 100 C	12
b. OST 100	12
c. OST 150	13
d. OST 250 C	13
e. OST 250	14
f. OST 250 H	14
g. OST 350	15
h. Standard features	15
2. List of accessories	16
a. Standard	16
b. Optional	16
c. Detachable parts	17
V. OPERATING INFORMATION	18
1. Installation of the device	19
2. Turning [On/Off] the device	19
3. Connection to other instruments	19
VI. USE OF THE DEVICE	21
1. OST 100C	22
2. OST 100, OST 150	22
3. OST 250 C	23
4. OST 250, OST 250 H, OST 350	23
5. Reminder & Best practices	24
a. OST012 - Manual phoropter arm	24
b. Motorized phoropter arm	27
c. Wheelchair base plate	27
d. Patient chair	28
e. Drawer	29
f. Power button	29

VII. MAINTENANCE	30
1. Storage and handling condition	31
2. Cleaning instructions	31
VIII. ERROR AND TROUBLESHOOTING	32
IX. TECHNICAL DESCRIPTION	34
1. Technical data	35
a. Product lifespan	35
b. Product dimensions and weight	35
c. LED specifications	35
d. Acoustic levels	35
e. Product electrical data	36
f. Power consumption	36
g. IP Classification	36
h. Technical Information	36
2. Electromagnetic compatibility	37
3. Disposal	39
X. SYMBOL EXPLANATION	40
1. On the document	41
2. On the device	41
3. On the packaging	42
XI. EXCLUSION OF LIABILITY	43
XII. QR CODE	45
XIII. CONTACT INFORMATION	49

I. INTRODUCTION





The latest version of this user manual is available on a web space.

To access other available languages, please scan the QR code available at the end of this user manual > QR Code Chapter (p.45).

For a safer, more effective use, follow the instructions outlined in this manual.

Copyright © 2026 Essilor - Original manual - All rights reserved.


Essilor International

147 rue de Paris, 94220, CHARENTON-LE-PONT

www.essilor.com

All reproduction of the content of this document, whether in part or as a whole, for the purpose of its publication or dissemination by any means and in any format whatsoever, even free of charge, is strictly prohibited without Essilor's prior written consent.

II. DIRECTION FOR USE



1. Intended use

a. Intended purpose

The OST Refraction Units are pieces of furniture intended to be used by eye care professionals to hold ophthalmic instruments (such as phoropter, refractometer, keratometer, slit lamp) and provide several patient and instruments ergonomical adjustments during ophthalmic examination.

2. Expected clinical benefit

These refraction units provide support and position control to ophthalmic instruments by eye care professionals and patients during ophthalmic examination or treatment.

3. Contraindications

There are no known contraindications.

These units are designed to accommodate patients up to 150 kg.

4. Side effects

There are no known side-effects.

Please report any serious incident that occurred in relation to the device to essilor-instruments-vigilance@essilor.com and to the local competent authority for medical devices.

5. Intended population


These refraction units are adapted to patients that can be seated on the chair, according to the units' maximum patient weight capacity.

The refraction units OST 250 H and OST 350 are designed to accommodate wheelchair patients.

6. Intended users

These devices are intended for eye care professionals use only.

III. CAUTIONS AND WARNINGS



1. Definitions

SYMBOL	DESCRIPTION
	Caution: a hazardous situation that, if not avoided, could result in minor or moderate injury.
	Warning: a hazardous situation that, if not avoided, could result in death or serious injury.
	Important and/or useful additional information to learn relating to the text in this manual.

2. Product safety

- Please do not use the refraction unit other than explained in this manual.
- Please use the refraction unit at 110-120 VAC or 220-240 VAC for OST100, OST100C, OST150, OST250C or 110-240 VAC 50-60 Hz for OST250, OST250H, OST350 (the tolerance range for these values has been determined as 10%). Otherwise, severe damage to the unit and the user may occur.

a. Precautions for use



- The covers are fragile, handling them while wearing jewellery or having long nails can lead to scratches.
- The white covers may yellow over time when exposed to ultraviolet light for an extended period.



- Never attempt to dismantle the instrument. This may cause a malfunction or fire.
- If the instrument does not work properly, do not touch the inside. Disconnect the plug from the outlet and consult your dealer.
- If liquid spills onto the instrument or foreign objects get inside, disconnect the plug from the outlet and consult your dealer.
- If any abnormalities occur (noise, smoke, etc.), disconnect the plug from the outlet and consult your dealer. Continued use may result in fire or personal injury.
- Do not set up the refraction unit in places where it may be leaking water or on a wet ground.
- After the installation, do not try to move the instrument table. Do not pull the console and column stand, otherwise damage to the unit may occur.
- Technical installation and set-up of the instrument should be made by authorized service personnel. Any interventions made to the refraction unit other than made by authorized service personnel shall put the unit out of warranty.
- During examinations, the unit's movements must be controlled under constant supervision by the healthcare professional. The healthcare professional must ensure that there are no obstructions during movement that could cause collisions leading to injury or equipment damage. The patient should be warned by the healthcare professional that the patient's hands should not come into contact with the table when the table is moved to the zero position.
- The healthcare professional must ensure that the space between the patient in the chair and the table of the unit to be moved is wide enough to prevent contact between the patient and the table. The healthcare professional must lower the patient chair and/or increase the table height before starting the examination.
- To handle the manual arm, always release the black handle **BEFORE** rotating the arm and then lock it to keep it in place.



- Do not try to repair or modify the instrument.
- Never try to perform any repairs inside the instrument yourself. In the event of malfunctions, consult your dealer.
- Children under the age of 6 should not use refraction unit without their parents.
- The refraction unit is not suited for use in air/flammable anesthetic gas, oxygen or nitrous oxide/flammable anesthetic gas atmosphere.
- This product should be used in an environment free of flammable anesthetic gas and other flammable gases.



- The patient should not get up from the chair or touch motorized arm while the motorized arm is opened. Otherwise, bumps, tripping, or accidents may occur.
- The patient should sit in or get out of the chair when it is in its lowest position, and the footrest is folded down. Otherwise, it may cause tripping, falls, and accidents.
- The base plate must be used after the refraction unit is fixed in place. The base plate must absolutely not be used before the refraction unit is fixed in place or after it has been removed from the refraction unit. Otherwise, it cause accidents and injuries.
- The base plate must absolutely not be used for any purpose after being removed from the refraction unit. Do not attempt to use or recline the base plate after it has been removed from the refraction unit. Otherwise, it may cause accidents and injuries.
- The patient should not get up from the chair while the drawer is open. Otherwise, bumps, tripping, or accidents may occur.

b. Power source



- Do not use multi-socket power strips, adapters or extension cords to connect the instrument to the mains.
- Make sure the power cord is fully inserted into both the outlet and the instrument. Failure to insert it properly may result in a fire or electric shock.
- If the power cord becomes hot after using the instrument, check that it is not dirty. If it is not, replace the power cord with a new one. Continued use may cause malfunction or personal injury.
- Use the instrument with the appropriate supply voltage. Continued use with a supply voltage greater than the rated power may cause malfunction or fire.
- Use only the power cord provided with the device.
- To avoid the risk of electric shock, this device must only be connected to a supply mains with an earth connection.
- Take care to use the power cord grounding cable when connecting to the ground terminal.
- Do not damage the power cord (by bending it, pulling it or placing heavy objects on top of it, etc.). Do not modify it either. If the cord is damaged (loose contact, damaged sheath, etc.), replace it with a new cord. Continued use may result in an electric shock or fire.
- Do not touch the power plug with wet hands. This may cause an electric shock.



If you do not use the instrument for an extended period, disconnect the power cord from the outlet.


c. Electromagnetic compatibility

- This device generates, uses and can radiate radio frequency (RF) energy. If this equipment is not used as specified in this manual, it may cause electromagnetic interference.
- This device has been tested in accordance with the EN 60601-1-2 Standard for Medical Products and its compliance with acceptable limits has been determined. These limits indicate that the device provides acceptable protection against electromagnetic interference (EMC) when used as directed in the manual.
- This device has been designed and manufactured in accordance with the requirements of the EN 60601-1-2 standard. It may be affected by portable and mobile RF communications devices. It should therefore not be stored with other equipment.
- For more information on this device and EMC, see part IX.2 of this manual.

d. Others

- The metal chassis of the refraction unit is protectively grounded. However, the reliability of the grounding depends on the electrical supply of the place where the unit is used. The unit should be used by plugging into an earthed socket.
- Please unplug the main power cable when working place lamp and socket power inlets on refraction unit are being plugged on and off.
- Check the refraction unit and accessories carefully and be sure they have not been damaged during dispatch. About the refraction unit and/or accessories that are damaged during dispatch, please get in touch with your dealer.
- Prior to first use of the refraction unit, please carefully read and understand the user manual thoroughly.
- Please keep the user manual at a distance where the user can reach quickly.

IV. PRODUCT DESCRIPTION



1. Product plan with description


a. OST 100 C



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Main body:** Main body of the refraction unit.
4. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.


b. OST 100



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Manual phoropter arm:** Arm where the phoropter is installed.
4. **Main body:** Main body of the refraction unit.
5. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.


c. OST 150



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Manual phoropter arm:** Arm where the phoropter is installed.
4. **Main body:** Main body of the refraction unit.
5. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.


d. OST 250 C



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Main body:** Main body of the refraction unit.
4. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.

e. OST 250



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Main body:** Main body of the refraction unit.
4. **Drawer:** Drawer to put the trial lens set or accessories in.
5. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.

f. OST 250 H



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Main body:** Main body of the refraction unit.
4. **Drawer:** Drawer to store the trial lens set or accessories.
5. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.

g. OST 350



With:

1. **Chair:** Chair where the patient is seated during examination.
2. **Stand column:** Column where the head lamp, chart projector, manual phoropter arm and other instruments are installed.
3. **Triple table:** Triple table that rotates, which can carry three ophthalmology instruments together.
4. **Drawer:** Drawer to store the trial lens set or accessories in.
5. **Main body:** Main body of the refraction unit.
6. **Plug outlet:** Plug outlet for an extra instrument 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country.

h. Standard features

Standard features	OST 100 C	OST 100	OST 150	OST 250 C	OST 250	OST 250 H	OST 350
Fixed chair backrest	✓	✓	✓	✓	NA	NA	NA
Fully reclinable chair backrest	NA	NA	NA	NA	✓	✓	✓
360° chair rotation	✓	✓	✓	NA	✓	✓	✓
Height-Adjustable motorized chair	✓	✓	✓	✓	✓	✓	✓
Separated baseplate	✓	NA	NA	NA	NA	NA	NA
Suitable for wheelchair	NA	NA	NA	NA	NA	✓	✓
5 W LED head lamp	✓	✓	✓	✓	✓	✓	✓
Flexible near vision lamp	Optional	Optional	Optional	Optional	✓	✓	✓
Chair footrest	✓	Optional	✓	✓	✓	✓	✓
Chair armrest	Optional	Optional	✓	Optional	✓	✓	✓
Manual arm	Optional	✓	✓	Optional	Optional	Optional	Optional
Sliding table	NA	✓	✓	✓	✓	✓	Rotating
Mini keypad on sliding table	NA	NA	NA	NA	✓	✓	NA
Mini keypad on rotating table	NA	NA	NA	NA	NA	NA	✓
Height adjustable electrical table	NA	NA	NA	✓	✓	✓	✓
Chart projector plate on profile	Optional	✓	✓	Optional	✓	✓	Optional
Projector arm plate	Optional	Optional	Optional	Optional	Optional	Optional	✓
Column accessory holder	Optional	Optional	Optional	Optional	Optional	✓	✓

Mechanical brake	NA	✓	✓	NA	✓	NA	NA
Electromagnetic brake	NA	NA	NA	✓	NA	✓	✓
Foot pedal for chair lift	NA	NA	NA	Optional	Optional	Optional	✓
Common chinrest	NA	Optional	Optional	Optional	✓	✓	Optional
Single drawer	NA	NA	NA	NA	✓	✓	✓

2. List of accessories

a. Standard

- User manual

b. Optional


Optional accessories	References	OST 100 C	OST 100	OST 150	OST 250 C	OST 250	OST 250 H	OST 350
Unit body RAL3001 side cover set	OST010	NA	✓	✓	NA	✓	NA	NA
Manual phoropter arm	OST012	✓	✓	✓	✓	✓	✓	✓
Unit body RAL5005 side cover set	OST020	NA	✓	✓	NA	✓	NA	NA
Unit body RAL9016 side cover set	OST030	NA	✓	✓	NA	✓	NA	NA
Spacer	OST046	✓	✓	✓	✓	✓	✓	✓
OST column accessory holder	OST047	✓	✓	✓	✓	✓	✓	✓
Monitor support arm	OST048	✓	✓	✓	✓	✓	✓	✓
Chair footrest	OST109	✓	✓	✓	✓	✓	✓	✓
Chair armrest	OST110	✓	✓	✓	✓	✓	✓	✓
Flexible near vision lamp	OST127	✓	✓	✓	✓	NA	NA	NA
Projector arm plate	OST202	✓	✓	✓	✓	✓	✓	✓
Flexible near vision lamp	OST208	NA	NA	NA	NA	✓	✓	✓
Common chinrest	OST210	NA	✓	✓	✓	✓	✓	NA
Unit body RAL5005 side cover set	OST251	NA	NA	NA	NA	NA	✓	NA
Unit body RAL9016 side cover set	OST252	NA	NA	NA	NA	NA	✓	NA
Unit body RAL3001 side cover set	OST253	NA	NA	NA	NA	NA	✓	NA
Double column	OST254	NA	NA	NA	NA	✓	✓	NA
Unit body RAL5005 side cover set	OST255	NA	NA	NA	✓	NA	NA	NA
Unit body RAL9016 side cover set	OST256	NA	NA	NA	✓	NA	NA	NA
Unit body RAL3001 side cover set	OST257	NA	NA	NA	✓	NA	NA	NA
Chair foot pedal	OST371	NA	NA	NA	NA	✓	✓	✓
Double drawer with shelf	OST376	✓	✓	✓	✓	✓	✓	✓
Motorized phoropter arm	OST380	✓	NA	NA	NA	✓	✓	✓

Mini motorized phoropter arm	OST380C	NA	NA	NA	✓	NA	NA	NA
Triple drawer mobile	OST381	✓	✓	✓	✓	✓	✓	✓
Double drawer mobile	OST382	✓	✓	✓	✓	✓	✓	✓
Unit body RAL5005 side cover set	OST384	NA	NA	NA	NA	NA	NA	✓
Unit body RAL9016 side cover set	OST385	NA	NA	NA	NA	NA	NA	✓
Unit body RAL3001 side cover set	OST386	NA	NA	NA	NA	NA	NA	✓
Depth adjustable seat mechanism	OST387	✓	✓	✓	✓	✓	✓	✓
Phoropter stand	OSTKBDESK	✓	✓	✓	✓	✓	✓	✓

c. Detachable parts

- Power cable (2.5 m)

V. OPERATING INFORMATION








- This instrument must be installed by authorized service personnel. To install the instrument or to change its connection, please contact your Essilor dealer.
- Respect the precautions below:
 - Do not install the instrument in a location:
 - Where dust or dirt accumulates,
 - Directly exposed to sun light,
 - Oxygen rich,
 - Displaying extreme temperatures and humidity levels,
 - Likely to undergo strong oscillations or sudden shocks.
- Do not use the instrument with flammable anesthetics or in conjunction with flammable agents

1. Installation of the device



Essilor declines all responsibility in the event of incorrect installation of the product by anyone other than an authorized service personnel trained by Essilor.

2. Turning [On/Off] the device

OST 100 C	OST 250 C	OST 100 / OST 150 / OST250 / OST250 H / OST 350
		



Make sure that you plug in the device.

[On/Off] operation is performed by pressing the power button indicated in red above.

3. Connection to other instruments


The OST can be connected to accessories, such as:

- Manual phoropter arm
- Spacer
- Column accessory holder
- Monitor support arm
- Chair footrest
- Chair armrest
- Flexible near vision lamp
- Projector arm plate
- Common chinrest
- Double column
- Chair foot pedal
- Double drawer with shelf
- Motorized phoropter arm
- Mini motorized phoropter arm
- Triple drawer mobile
- Double drawer mobile

For the other instruments, technical installation and set-up must be made by authorized service personnel.

For details, please ask your local dealer.

VI. USE OF THE DEVICE



1. OST 100C

Keypad	Control	Description
		Power On-Off control
		Chair Up-Down control
		(Optional) motorized phoropter arm control
		Near vision lamp control

2. OST 100, OST 150

Keypad	Control	Description
		Power On-Off control
		Lamp light control
		Projector power control
		Auxiliary power control
		Chair Up-Down control
		Slit lamp light control

3. OST 250 C

Keypad	Control	Description
		Power On-Off control
		Lamp light control
		Projector power control
		Auxiliary power control
		Chair Up-Down control
		Slit lamp light control
		Electro magnet control
		Table Up-Down control
		(Optional) Motorized phoropter arm control

4. OST 250, OST 250 H, OST 350


Keypad	Control	Description
		Power On-Off control
		Dimmable head lamp control
		Projector power control
		Auxiliary power control
		(Optional) Motorized phoropter arm control
		Near vision lamp control
		Room lamp control
		Electro magnet control
		Chair Up-Down control

		Table Up-Down control
		Slit lamp light control

5. Reminder & Best practices

a. OST012 - Manual phoropter arm



To prevent the vertical movement and rotation

Use the black handle (c) to block vertical position and rotation.



The black handle (c) should be loosened before performing the movement and tightened again after the operation completed.




The rotation movement should not be performed while the black handle (c) is fully tightened.




> Phoropter arm is fixed.

Incorrect usage	Correct usage
Do not push the phoropter arm while the fixing knob (c) is fully tightened.	Loosen the fixing knob before moving the phoropter arm.
	




Adjust of vertical position (to lock & unlock)

Use the handle (h) for blocking or unblocking.




While moving the phoropter arm up and down safely and in a controlled manner, support the phoropter arm with your hand.




Adjust of horizontal direction (to lock & unlock)

Use the regulator (e) to adjust properly horizontal direction.




Use the regulator (d) to displace the phoropter in horizontal position.



To change the regulation force:


- 1 Loosen the safety screw (b).
- 2 With an Allen key, turn the screw (a) to get more or less force F.
- 3 Tighten the safety screw (b).



b. Motorized phoropter arm




The patient should not get up from the chair or touch motorized phoropter arm while the motorized arm is opened.
Otherwise, bumps, tripping, or accidents may occur.



c. Wheelchair base plate



The patient should sit in or get out of the chair when it is in its lowest position, and the footrest is folded down.
Otherwise, it may cause tripping, falls, and accidents.



d. Patient chair



The base plate must be used after the refraction unit is fixed in place.

The base plate must absolutely not be used before the refraction unit is fixed in place or after it has been removed from the refraction unit.

Otherwise, it may cause accidents and injuries.



The base plate must absolutely not be used for any purpose after being removed from the refraction unit. Do not attempt to use or recline the base plate after it has been removed from the refraction unit.

Otherwise, it may cause accidents and injuries.




When adjusting the seat back to recline position, support the backrest with your hand to ensure the seat reclines safely and in a controlled manner.





Be careful not to hit the refraction unit when turning the seat.



e. Drawer



The patient should not get up from the chair while the drawer is open.
Otherwise, bumps, tripping, or accidents may occur.



f. Power button




The refraction unit can be controlled during use with the power button on the keypad.

At the end of the day or if it is not used for a long time, the refraction unit should be switched off using the power button on the main body.

Otherwise, it will cause unnecessary energy consumption.



VII. MAINTENANCE





In order to ensure safety and the performance of the instrument, all maintenance operations, unless otherwise specified in this manual, must be carried out by qualified maintenance technicians.



If you find that this device is dirty, you can clean it as often as you want (see after the specific cleaning methods).

1. Storage and handling condition



Respect the operating, storage and transport conditions noted below.

Avoid condensation conditions.

	Temperature	Humidity	Atmospheric pressure
Use	[+20°C; +30°C]	[30 %; 70%]	[0.49 atm; 1 atm]
Storage	[0°C; +40°C]	[30 %; 70%]	[0.49 atm; 1 atm]
Transport	[-10°C; +50°C]	[30 %; 70%]	[0.49 atm; 1 atm]

2. Cleaning instructions




Never use alcohol-xylol and similar liquid cleaner for the cleaning of the refraction unit chair and wooden part. These cleaning materials may cause damage on the plastic and wooden surfaces.



To avoid any incident, unplug the instrument before cleaning.

- Prior to cleaning, please be sure to turn off the device by pressing the power button and to unplug the power cable.
- You can clean both the wooden and metal parts of the refraction unit with a clean damp cloth.

VIII. ERROR AND TROUBLESHOOTING



If a problem is detected, refer to the table below in order to take the appropriate measures.

Symptoms	Causes and Measurements
Refraction unit does not work. Power keypad led does not turn on.	No electricity on the unit - The outlet where the instrument is connected has a breakdown. Call an authorized electrician to check the outlet. - Power cable of the table is damaged (signs of being smashed or cracked can be seen). Spare power cable can be provided by your authorized dealer. - Check the glass fuse at the back of the instrument. Replace with fuses with same value. Spare fuses can be provided by your authorized dealer. Keypad panel does not work - Contact authorized service.
Head lamp does not work	Head lamp bulb has blown Turn off power, unplug the cable and change the bulb. Spare bulbs can be provided by your authorized dealer.


If the problem has not been resolved after taking the measures listed above, contact your local distributor immediately.
Your dealer has been trained by Essilor.



Please refer to the name plate and let us have the following information:

- Name of the instrument: OST XXX
- Serial number: 7 digit-characters indicated on the name plate
- Detailed defects

IX. TECHNICAL DESCRIPTION



OST is a class I medical device (Regulation (EU) 2017/745).

Basic UDI of the device: 361502000000I0ST0000000DX.

1. Technical data

a. Product lifespan

The expected life of the device and its components is 10 years.

b. Product dimensions and weight

Dimensions

Dimensions	OST 100 C	OST 100	OST 150	OST 250 C	OST 250	OST 250 H	OST 350
Height	195 cm	185 cm	185 cm	185 cm	185 cm	185 cm	180 cm
Width	Body: 50 cm Chair: 50 cm	160 cm	160 cm	140 cm	175 cm	180 cm	175 cm
Depth	Body: 67.5 cm Chair: 89 cm	130 cm	130 cm	118.5 cm	142 cm	160 cm	185 cm

Weight

Weight	OST 100 C	OST 100	OST 150	OST 250 C	OST 250	OST 250 H	OST 350
Total weight	160 kg	190 kg	200 kg	245 kg	235 kg	375 kg	360 kg

c. LED specifications

For OST 100 C, OST 100, OST 150, OST 250 C

- 3-5 W
- Warm white
- 110-120/220-240 VAC, 50-60 Hz
- Type: GU10

For OST 250, OST 250 H, OST 350

- 3-5 W Dimmable
- Warm white
- 110-120/220-240 VAC, 50-60 Hz
- Type: GU10

d. Acoustic levels

For OST 100 C, OST 100, OST 150, OST 250 C

- Ambient: 45 dB
- Maximum exposure: 53 dB

For OST 250, OST 250 H, OST 350

- Ambient: 48.2 dB
- Maximum exposure: 58.2 dB

e. Product electrical data

For OST 100C

- Power: 110-120 VAC or 220-240 VAC, 50-60 Hz
- Fuse: 8 Amper
- Electrical classification: Class 1, Type B

For OST 100, OST 150, OST 250 C

- Power: 110-120 VAC or 220-240 VAC, 50-60 Hz
- Fuse: 8 Amper
- Table electricity outlet: 0-6-12 VAC Dimmable and 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country
- Electrical classification: Class 1, Type B

For OST 250, OST 250 H

- Power: 110-240 VAC, 50-60 Hz
- Fuse: 8 Amper
- Table electricity outlet: 0-6-12 VAC Dimmable and 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country
- Electrical classification: Class 1, Type B

For OST 350

- Power: 110-240 VAC, 50-60 Hz
- Fuse: 10 Amper
- Table electricity outlet: 0-6-12 VAC Dimmable and 110-120 VAC or 220-240 VAC, 50-60 Hz depending on customer country
- Electrical classification: Class 1, Type B

f. Power consumption

Model	Volt	Amper	Volt Amper	Watt	Kilowatt
OST 100 C	110 V	1,4 A	110 VA	100 W	1 KWH
OST 100/OST 150	110 V	1,4 A	110 VA	100 W	1 KWH
OST 250 C	110 V	2,4 A	210 VA	190 W	1,8 KWH
OST 250	110 V	1,2 A	102 VA	85 W	1 KWH
OST 250 H	110 V	1,4 A	110 VA	100 W	1 KWH
OST 350	110 V	1,8 A	150 VA	140 W	1,3 KWH

g. IP Classification

IP20

h. Technical Information

	OST 100 C	OST 100	OST 150	OST 250 C	OST 250	OST 250 H	OST 350
Sliding table Capacity	NA	20+20 kg	20+20 kg	20+20 kg	20+20 kg	20+20 kg	20+20+20 kg
Chair load capacity	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg	150 kg
Chair fixed by hand	360° rotation	360° rotation	360° rotation	Fixed	360° rotation	360° rotation	360° rotation
Chair stroke	18 cm	18 cm	18 cm	18 cm	18 cm	18 cm	18 cm
Chair motor	24 V DC	24 V DC	24 V DC	24 V DC	24 V DC	24 V DC	24 V DC
Table stroke	NA	NA	NA	18 cm	18 cm	18 cm	18 cm
Table motor	NA	NA	NA	24 V DC	24 V DC	24 V DC	24 V DC

2. Electromagnetic compatibility

OST Refraction Units conforms to the requirements of the EMC (electromagnetic compatibility) standard.

This device complies with EMC standard IEC 60601-1-2: 2014+A1:2020, and the expected electromagnetic environment for the entire life cycle is the Occupational healthcare facility environment.



- The device complies with the applicable electromagnetic compatibility standards, however, the user must ensure that any electromagnetic interference does not create an additional risk, such as radio frequency transmitters or other electronic devices.
- The different device's cords must be separated from each other.
- Certain types of mobile telecommunication devices, such as mobile phones, may interfere with the device. Recommended separation distances must therefore be respected.
- The device shall not be used in the vicinity of or placed on another device. If this cannot be avoided, it is necessary to check its proper functioning under the conditions of use before using it. The use of accessories other than those specified or sold by the manufacturer as replacement parts may result in an emissions increase or a decrease in the immunity of the device.
- In case the device stops working, reset the device and restart test from the beginning. Do not use the previous data to make prescription.
- When using this device in hospitals, do not place it near active HF surgical equipment or in RF shielded rooms with an ME system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.

Guidance and manufacturer's declaration - electromagnetic emissions

The OST Refraction Unit is intended for use in the electromagnetic environment specified below. The customer or user of the OST Refraction Unit should assure that it is used in such an environment.

Emissions test	Basic EMD standard	Compliance
Conducted and radiated RF emissions	CISPR 11	Class B, Group 1
Harmonic distortion	IEC 61000-3-2	Class A
Voltage fluctuations and flicker	IEC 61000-3-3	Complies

Immunity test	Basic EMC standard or test method	Immunity test levels Home healthcare environment	Compliance level
Electrostatic discharge	IEC 61000-4-2	± 8kV contact ± 2, 4, 8, 15kV air	± 8kV contact ± 2, 4, 8, 15kV air
Radiated RF EM fields	IEC 61000-4-3	3 V/m ^a 80 MHz to 2.7 GHz 80% MA at 1 kHz	3V/m
Proximity fields from RF wireless communications equipment		See the table below.	
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m
Proximity magnetic fields	IEC 61000-4-39	30kHz (8A/m) 134.2kHz (65A/m) 13.56 MHz (7.5A/m)	30kHz (8A/m) 134.2kHz (65A/m) 13.56MHz (7.5A/m)

^a Before modulation is applied.

Test frequency (MHz)	Band ^a (MHz)	Service ^a	Modulation ^b	Maximum power (W)	Distance (m)	Immunity test level (V/m)	Compliance level (V/m)
385	380 - 390	TETRA400	Pulse modulation ^b 18Hz	1.8	0.3	27	27
450	430 - 470	GMRS460, FRS460	FM ±5kHz deviation 1kHz sine	2	0.3	28	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^b 217Hz	0.2	0.3	9	9
745							
780							
810	800 - 960	GSM800/900, TETRA800, iDEN820, CDMA850, LTE Band 5	Pulse modulation ^b 18Hz	2	0.3	28	28
870							
930							
1720	1700 - 1990	GSM1800 ; CDMA1900 ; GSM1900 ; DECT ; LTE Band 1, 3, 4, 25 ; UMTS	Pulse modulation ^b 217Hz	2	0.3	28	28
1845							
1970							
2450	2400 - 2570	Bluetooth, WLAN, 802.11b/g/n, FRID2450, LTE Band 7	Pulse modulation ^b 217Hz	2	0.3	28	28
5240	5100 - 5800	WLAN 802.11a/n	Pulse modulation ^b 217Hz	0.2	0.3	9	9
5500							
5785							

^a For some services, only the uplink frequencies are included.

^b The carrier shall be modulated using a 50% duty cycle square wave signal.

Immunity test	Basic EMC standard	Immunity test levels Home healthcare environment	Compliance level
Electrical fast transients/ bursts	IEC 61000-4-4	Input AC power port ± 2kV 100 kHz repetition frequency	± 2kV
		Signal input/ output unit port ±1kV 100 kHz repetition frequency	± 1kV
Surges Line-to-line	IEC 61000-4-5	± 0.5kV, 1kV	
Surges Line-to-ground		± 0.5kV, ± 1kV, ± 2kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 Vrms 0.15 MHz - 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	3 Vrms 6Vrms in ISM bands

Voltage dips	IEC 61000-4-11	0% U_T ; 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% U_T ; 0.5 cycle
		0% U_T ; 1 cycle and 70% U_T ; 25 cycle Single phase: 0°	0% U_T ; 1 cycle 70% U_T ; 25 cycle
Voltage interruptions		0% U_T ; 250 cycle	0% U_T ; 250 cycle

U_T is the AC supply voltage before applying inspection level.


3. Disposal

	Instructions for the disposal of the instrument in accordance with Directives 2012/19/EU and 2011/65/EU regarding the limitation of dangerous substances in electrical and electronic equipment and the disposal of electrical and electronic waste. When it reaches the end of its lifetime, the instrument should not be thrown out with the household refuse. It can be disposed of at a waste management center operated by the municipality or the retailers who offer this service. The separate disposal of an electrical device avoids any damage to the environment or health that could result from a non-compliant disposal, and also allows the materials it is composed of to be recycled in order to save energy and resources. The pictogram of the wheeled container appears on the label of the instrument. It indicates the obligation for separate collection and disposal of end-of-life/out-of-use electrical and electronic equipment.
---	--



- The user must take into account the potentially harmful effects on the environment and human health that could result from the non-compliant disposal of the instrument in its entirety or some of its components.
- To avoid the release of dangerous substances into the environment and to encourage the preservation of natural resources, the manufacturer facilitates, in the event that the user wishes to dispose of the instrument at the end of its lifespan, the reuse, recovery and recycling of the instrument and its components. Before disposing of the instrument, the requirements of European and national regulations must be taken into consideration.
- Do not dispose of the instrument with household waste, but dispose of it separately by giving it in a company specialized in the disposal of electrical and electronic equipment or at the local administrative services in charge of waste collection.
- The supplier or manufacturer is required to recover the old equipment.
- By joining a consortium for the waste of technological equipment, the manufacturer covers the treatment and recycling costs of the used instrument.
- The manufacturer undertakes to provide the user with all the information relating to the dangerous substances contained in the device and the methods of recycling these substances, and to inform them of the existence of recycling of the used equipment. The law provides for severe penalties in case of infringement.

X. SYMBOL EXPLANATION



1. On the document

SYMBOL	DESCRIPTION
	Caution: a hazardous situation that, if not avoided, could result in minor or moderate injury.
	Warning: a hazardous situation that, if not avoided, could result in death or serious injury.
	Important and/or useful additional information to learn relating to the text in this manual.
	Tips: practical advice.

2. On the device


SYMBOL	DESCRIPTION
	Manufacturer
	Country of manufacture and Manufacturing date (Year-Month)
	Medical device
	CE Marking (European regulation relating to medical devices)
	Product reference
	Product serial number
	Unique device identifier
	Refer to instruction manual/ booklet
	General warning sign
	Alternating current
	Protective earth (ground)
	Type B applied part
	IP classification
	Waste disposal symbol in accordance with Directives 2012/19/EU and 2011/65/EU

3. On the packaging

For proper handling, storage and transport requirements.

SYMBOL	DESCRIPTION
	This way up
	Maximum stacking of products above market product
	Fragile; handle with care
	Keep dry
	Indicate the thermal limits to which the medical device can be exposed in complete safety
	Indicate the humidity limits to which the medical device can be exposed in complete safety
	Indicate the limits of atmospheric pressure to which the medical device can be exposed in complete safety
	Medical device
	CE Marking (European regulation relating to medical devices)
	Product reference
	Product serial number
	Unique device identifier
	Waste disposal symbol in accordance with Directives 2012/19/EU and 2011/65/EU
	Indicates that the packaging is recyclable

XI. EXCLUSION OF LIABILITY



The product shall be used in compliance with the applicable laws and regulations, by qualified, professional users. The product must be installed and used according to the instructions provided in the present user manual and to any written direction or recommendation provided by Essilor (the “documentation”).

Essilor reserves the right to revise the documentation and to make changes in its content from time to time. Preventive and corrective maintenance (including regular calibration if necessary, according to the documentation) shall be performed in accordance with the documentation.


Any product warranty offered by Essilor is contingent upon use of the product in accordance with the documentation and with the product's intended use and does not cover products that were modified without Essilor's prior written approval or repaired by a third party not approved by Essilor, nor to products that were subjected to physical, chemical or electrical stress that the products were not originally designed for.

Essilor shall not be held liable for any damages suffered by the user of the product, the product, or any third party, resulting from any user's non-compliance to the present section.

If the product offers a connectivity function, the user shall be sole responsible for:

- selecting, obtaining and maintaining all required internet access and telecommunications at its own costs; and
- adopting and maintaining procedures and measures to protect its workstations, hardware and software, other than the Product, including against any virus or intrusion

XII. QR CODE




The latest version of the user manual in the appropriate language is available on a web space. Upon request, a paper version can be provided for free.

en	The complete user manual is available on a web space in PDF format. To access it, please scan the QR code below using a dedicated tool or application. Please make sure that your device is suitable and has an appropriate software to display the electronic Instructions for use.
fr	Le manuel utilisateur complet est disponible sur un espace web au format PDF. Pour y accéder, veuillez scanner le QR code ci-dessous à l'aide d'un outil ou d'une application dédié(e). Veuillez vous assurer que votre appareil est compatible et dispose d'un logiciel approprié pour afficher le manuel électronique.
ar	لتمكن من الوصول إليه، يُرجى مسح رمز الاستجابة السريعة PDF. دليل المستخدم الكامل متوفر من خلال موقع الويب بصيغة أدناه باستخدام أداة أو تطبيق مخصص لذلك. يُرجى التأكد من أن جهازك مناسب ويحتوي على برنامج مناسب لعرض التعليمات الإلكترونية الخاصة بالاستخدام.
be	Поўная інструкцыя карыстальніка даступна ў інтэрнэт-прасторы у фармаце PDF. Каб атрымаць да яе доступ, адсканіруйце QR-код ніжэй пры дапамозе спецыяльнага сродку або праграмы. Калі ласка, упэўніцеся, што ваша прылада прыдатная для паказу электроннай Інструкцыі па карыстанню і што на ёй усталявана адпаведнае праграмае забеспячэнне.
bg	Пълното ръководство за потребителя е достъпно в уеб пространството. За да получите достъп до него, моля, сканирайте QR кода по-долу, като използвате специален инструмент или приложение. Моля, уверете се, че вашето устройство е подходящо и разполага с подходящ софтуер за преглед на електронните Инструкции за употреба.
cs	Kompletní uživatelský návod je k dispozici na webovém prostoru ve formátu PDF. Chcete-li k němu získat přístup, naskenujte prosím níže uvedený QR kód pomocí speciálního nástroje nebo aplikace. Ujistěte se prosím, že používáte vhodné zařízení, které má vhodný software pro zobrazení elektronického uživatelského návodu.
da	Den komplette brugervejledning er tilgængelig på et webområde i PDF-format. For at få adgang til den skal du scanne QR-koden nedenfor ved hjælp af et dedikeret værktøj eller program. Sørg for, at din enhed er egnet og har en passende software til at vise de elektroniske brugsanvisninger.
de	Die vollständige Bedienungsanleitung ist auf einem Webspace im PDF-Format verfügbar. Für den Zugriff scannen Sie bitte den untenstehenden QR-Code mit einem speziellen Tool oder einer Anwendung. Bitte vergewissern Sie sich, dass Ihr Gerät für die Anzeige der elektronischen Gebrauchsanweisungen geeignet ist und über eine entsprechende Software verfügt.
el	Το πλήρες εγχειρίδιο χρήσης είναι διαθέσιμο σε έναν ιστοχώρο σε μορφή PDF. Για να αποκτήσετε πρόσβαση σε αυτό, σκανάρετε τον κωδικό QR παρακάτω χρησιμοποιώντας ένα ειδικό εργαλείο ή εφαρμογή. Βεβαιωθείτε ότι η συσκευή σας είναι κατάλληλη και έχει το κατάλληλο λογισμικό για την προβολή των ηλεκτρονικών οδηγιών χρήσης.
es	El manual de uso completo está disponible en un espacio web. en formato PDF. Para acceder a él, escanee el código QR debajo utilizando una herramienta o aplicación dedicada. Asegúrese de que su dispositivo sea adecuado y tenga el software apropiado para mostrar las Instrucciones de uso electrónicas.
et	Täielik kasutusjuhend on saadaval veebis PDF-vormingus. Juurdepääsuks palun skannige allolevat QR-koodi, kasutades selleks vastavat tööriista või rakendust. Veenduge, et teie seade sobib ja et selles on elektroonilise kasutusjuhendi kuvamiseks sobiv tarkvara.
fi	Täysi käyttöopas on saatavana verkosta PDF-muodossa. Saat pääsyn siihen skannaamalla alla olevan QR-koodin käyttäen siihen tarkoitettu työkalua tai sovellusta. Varmista, että laitteesi on sopiva ja sisältää asianmukaisen ohjelmiston sähköisten käyttöohjeiden esittämiseen.
he	למטה באמצעות כלי או QR-כדי לגשת אליו, יש לסרוק את קוד ה-PDF המדריך המלא למשתמש זמין באתר אינטרנט בפורמט אפליקציה ייעודיים. חשוב לוודא שהמכשיר שלך מתאים ובעל תוכנה מתאימה להצגת הוראות השימוש האלקטרוניות.
hr	Potpun korisnički priručnik dostupan je na mrežnom prostoru u PDF formatu. Da biste mu pristupili, skenirajte QR kod u nastavku pomoću odgovarajućeg alata ili aplikacije. Provjerite je li vaš uređaj prikladan i ima li odgovarajući softver za prikaz elektroničkih uputa za upotrebu.


hu	A teljes felhasználói kézikönyv elérhető az interneten PDF formátumban. Eléréséhez olvassa be az alábbi QR-kódot egy erre szolgáló eszközzel vagy alkalmazással. Ellenőrizze, hogy eszköze képes és rendelkezik a megfelelő szoftverrel az elektronikus használati útmutató megjelenítésére.
id	Panduan pengguna lengkap tersedia di ruang web dalam format PDF. Untuk mengaksesnya, silakan pindai kode QR di bawah ini menggunakan alat atau aplikasi khusus. Pastikan peranti Anda sesuai dan memiliki perangkat lunak yang layak untuk menampilkan petunjuk penggunaan elektronik.
it	Il manuale utente completo è disponibile in formato PDF su uno spazio Web. Per accedervi, leggere il codice QR sottostante mediante un apposito strumento o un'applicazione dedicata. Assicurarsi che il dispositivo sia adatto e che disponga di un software appropriato per visualizzare le istruzioni per l'uso in formato elettronico.
ja	完全なユーザーマニュアルは、PDF形式でウェブスペースから入手できます。アクセスするには、専用のツールまたはアプリケーションを使用して、以下のQRコードをスキャンしてください。お使いのデバイスが適切であり、電子説明書を表示する適切なソフトウェアがインストールされていることを確認してください。
ko	전체 사용 설명서는 웹 공간에 PDF 형식으로 있습니다. 이 설명서에 액세스하려면, 전용 도구 또는 앱을 사용하여 아래 QR 코드를 스캔하십시오. 사용자의 기기가 적합하고 전자적인 사용 설명서를 표시할 수 있는 적절한 소프트웨어가 있는지 확인하시기 바랍니다.
lt	Išsamaus naudotojo vadovo PDF formatu ieškokite interneto svetainėje. Kad jį atvertumėte, specialiu įrankiu arba programėle nuskaitykite toliau pateiktą QR kodą. Įsitikinkite, kad jūsų įrenginys yra tinkamas ir turi tinkamą programinę įrangą elektroninems naudojimo instrukcijoms rodyti.
lv	Pilnā lietotāja instrukcija ir pieejama tīmeklī PDF formātā. Lai tai piekļūtu, lūdzu, noskenējiet tālāk redzamo kvadrātkodu, izmantojot tam paredzētu rīku vai lietojumprogrammu. Lūdzu, pārlicinieties, vai jūsu ierīce ir piemērota un vai tai ir atbilstoša programmatūra elektroniskās lietotāja instrukcijas attēlošanai.
ms	Manual pengguna yang lengkap boleh didapati di ruang laman dalam format PDF. Untuk mengaksesnya, sila imbas kod QR di bawah menggunakan alat atau aplikasi khusus. Sila pastikan yang peranti anda adalah serasi dan mempunyai perisian yang sesuai untuk memaparkan Arahan elektronik untuk tujuan penggunaan.
mt	Il-manwal tal-utent shiĥ huwa disponibbli fuq il-web f'format PDF. Biex taċċessah, jekk jogħġbok skennja l-kodiċi QR t'hawn taht permezz ta' għodda jew applikazzjoni apposta. Jekk jogħġbok żgura li l-apparat huwa xieraq u għandu s-software adattat biex juri l-Istruzzjonijiet għall-Użu elettronici.
nl	De volledige gebruikershandleiding is in PDF-formaat beschikbaar op een website. U kunt de handleiding bereiken door de QR-code hiernaast te scannen met een geschikte applicatie. Uw apparaat moet geschikt zijn en over de juiste software beschikken om de elektronische gebruiksaanwijzing weer te geven.
no	Den komplette brukerhåndboken er tilgjengelig på et webhotell i PDF-format. For å få tilgang til den, skann QR-koden nedenfor ved hjelp av et dedikert verktøy eller applikasjon. Sørg for at enheten din er egnet og har en passende programvare for å vise den elektroniske bruksanvisningen.
pl	Kompletna instrukcja użytkownika jest dostępna na stronie internetowej w formacie PDF. Aby uzyskać dostęp, zeskanuj poniższy kod QR przy użyciu dedykowanego narzędzia lub aplikacji. Upewnij się, że urządzenie jest zgodne i wyposażone w odpowiednie oprogramowanie pozwalające wyświetlać elektroniczną Instrukcję obsługi.
pt	O manual do utilizador completo está disponível num espaço online no formato PDF. Para aceder a este, queira digitalizar o QR Code abaixo usando uma ferramenta ou uma aplicação dedicada. Certifique-se de que o seu dispositivo é compatível e possui um software apropriado para exibir as instruções eletrônicas de utilização.
pt (brazil)	O manual do usuário completo está disponível em um espaço online no formato PDF. Para acessar a este, por favor, digitalizar o QR Code abaixo usando uma ferramenta ou um aplicativo dedicado. Seu dispositivo deve ser compatível e possuir um software apropriado para exibir as instruções eletrônicas de utilização.
ro	Manualul de utilizare complet este disponibil online în format PDF. Pentru a-l accesa, scanați codul QR de mai jos folosind un instrument sau o aplicație dedicată. Asigurați-vă că dispozitivul dumneavoastră este potrivit și are un software adecvat pentru afișarea Instrucțiunilor de utilizare în format electronic.
ru	Полное руководство пользователя доступно в интернет-пространстве в формате PDF. Чтобы получить к нему доступ, отсканируйте QR-код ниже с помощью специального инструмента или приложения. Убедитесь, что

ваше устройство подходит и имеет соответствующее программное обеспечение для отображения электронных инструкций по эксплуатации.

- sk Celý používateľský manuál je dostupný vo webovom priestore vo formáte PDF. Ak chcete získať prístup, naskenujte nižšie uvedený QR kód pomocou špeciálneho nástroja alebo aplikácie. Uistite sa, že máte vhodné zariadenie s vhodným softvérom na zobrazenie elektronického návodu na použitie.
- sl Celoten uporabniški priročnik je na voljo kot dokument PDF na spletnem mestu. Za dostop optično preberite spodnjo kodo QR z namenskim orodjem ali aplikacijo. Prepričajte se, da je vaša naprava primerna in ima ustrezno programsko opremo za prikaz elektronskih navodil za uporabo.
- sr Kompletно uputstvo za korisnike je dostupno na veb prostoru u PDF formatu. Da biste mu pristupili, skenirajte QR kôd u nastavku pomoću namenske alatke ili aplikacije. Proverite da je vaš uređaj odgovarajući i da li ima potreban softver za prikaz elektronskog Uputstva za upotrebu.
- sv Den fullständiga bruksanvisningen finns tillgänglig på ett webbutrymme i PDF-format. För att komma åt den, vänligen skanna QR-koden nedan med ett dedikerat verktyg eller program. Se till att din enhet är lämplig och har en passande programvara för att visa de elektroniska användningsinstruktionerna.
- th สามารถรับคู่มือผู้ใช้ฉบับสมบูรณ์ในรูปแบบ PDF ได้จากบนเว็บไซต์ โดยในการเข้าถึง โปรดสแกนคิวอาร์โค้ดด้านล่างด้วยเครื่องมือหรือแอปพลิเคชันเฉพาะ โปรดตรวจสอบให้แน่ใจว่าอุปกรณ์ของคุณนั้นเหมาะสม และมีซอฟต์แวร์ที่สามารถใช้ในการแสดงคำแนะนำการใช้งานอิเล็กทรอนิกส์ได้อย่างถูกต้อง
- tr Kullanım kılavuzunun tamamı web alanında, PDF formatında mevcuttur. Buna erişmek için lütfen uygun bir araç veya uygulama kullanarak aşağıdaki QR kodunu okutun. Lütfen cihazınızın uyumlu ve elektronik kullanım talimatlarını görüntülemek için uygun bir yazılıma sahip olduğundan emin olun.
- uk Повна версія посібника користувача доступна в інтернеті в форматі PDF. Щоб отримати до нього доступ, скануйте QR-код нижче за допомогою спеціального додатку. Для перегляду електронного посібника користувача на вашому пристрої він повинен мати відповідні характеристики та програмне забезпечення.
- vi Hướng dẫn sử dụng đầy đủ có sẵn trên không gian web ở định dạng PDF. Để truy cập, vui lòng quét mã QR bên dưới bằng công cụ chuyên dụng hoặc bằng ứng dụng. Vui lòng đảm bảo rằng thiết bị của bạn phù hợp và có phần mềm phù hợp để hiển thị Hướng dẫn sử dụng điện tử
- zh 完整的操作手册以 PDF 格式在网络上提供。如需获取，请使用专门的工具或应用程序扫描下方二维码 QR 条码。请确保您的设备适用并安装有相应的软件，能够显示电子版使用说明。



XIII. CONTACT INFORMATION



If the instrument appears malfunctioning, it is highly recommended to check the instrument according to the troubleshooting procedure of this manual.

If any problem persists or the instrument is damaged or malfunctioning or it is mentioned to contact your local distributor, please follow the steps below.

- Please contact the local distributor in your province or country at first. All information is available on www.essilor-instruments.com in the "Contact" section.
- If the product has been provided with electronic instruction and you need a paper format, please contact your local distributor.
- If any serious incident occurred in relation to the device, report it to essilor-instruments-vigilance@essilor.com and to the local competent authority for medical devices.
- Before calling to the local distributor, please be sure to check Model and Serial Numbers.
- The serial number is unique to this unit and is accessible on the product. It is recommended to fill up the following table as soon as you purchase our product.
- Please keep this manual as a permanent record of your purchase and keep your purchase receipt as your proof of purchase.

Date of Purchase:

Dealer's Name:

Dealer address:

Dealer Phone No.:

Model No.:

Serial No.:



Essilor International
147, rue de Paris – 94220 Charenton-le-Pont France
www.essilor.com

